

Sample stock

meniscus polished

Artikel-Nr. / Item-No.	Bezeichnung / type	Ø mm	e / efl mm	R mm	Md. / ct. mm	material	Besonderheiten / specials	Menge / quantity
I.0100.0401	meniscus	4		6,61/-7,7	4	N- LASF31	AR @ vis area	2
I.0101.0401	meniscus	4		6,61/-7,7	4	N- LASF31	AR @ vis area	6
I.0110.0601	meniscus	6	38,21	16,600/- 100,930	1,5	BK7		10
I.0110.0802	meniscus	8		9857/- 12,47	3	F2		2
I.0111.0904	meniscus	9		7,5/-5,0	1,2	Borofloat	coated	1
I.0111.0002	meniscus	11		-0,563	0,5	N-LAK8		3
xxxxxxxxx	meniscus	12	34,5	15,90/- 46,42		SF10		2
I.0110.1203	meniscus	12	31,75	13,050/- 60,120	2,2	BK7		4
I.0111.1202	meniscus	12		R = -35,47 / 11,64	2,25	NSK2	AR @ 1060 nm	11
I.0111.1201	meniscus	12	52,82	16,300/- 37,920	3	BK7		0
I.0110.1301	meniscus	13		31,65/-20,7	1	SF6		1
I.0111.1301	meniscus	13	-97,26	20,680/- 14,000	2,7	BK7		0
I.0111.1411	meniscus	14,1		17,719/- 7,608	1,51	SK16		5
I.0101.1701	meniscus	17		20,419/- 37,925	2,5	LaK9	AR @ 904nm	1
I.0101.1741	meniscus	17,4		14,2/-17,5		SF14	AR @ 670 nm	4
I.0110.1781	meniscus	17,8		31,64/-20		SF6		3
I.0110.1801	meniscus	18		9,253/-118	9,6	SF2		4
I.0111.1802	meniscus	18		10,4/-95,09	6	SF2	coated	1
I.0111.1803	meniscus	18		126,70/- 130,0	2	B270	AR @ 700-1000 nm	3
I.0111.1804	meniscus	18		55,431/- 20,983	3	SF6	coated	2
I.0101.1901	meniscus	19		23,9/-4,2	3	SF14	AR @ 670 nm	2
I.0111.1902	meniscus	19		42,07/- 88,04	3,99	BK7		1
I.0101.2001	meniscus	20		20,67/- 33,48	4,5	SF11	AR @ 800-940 nm	1
I.0110.2243	meniscus	22,4	-46,61	37,480/- 12,470	17	B270		1



Sample stock

meniscus polished

Artikel-Nr. / Item-No.	Bezeichnung / type	Ø mm	e / efl mm	R mm	Md. / ct. mm	material	Besonderheiten / specials	Menge / quantity
I.0110.2244	meniscus	22,4	-43,96	34,750/- 12,470	10	B270		0
I.0110.2245	meniscus	22,4		36,54/-12,0	17	B270		10
I.0110.2247	meniscus	22,4	-47,06	37,480/- 12,470	17	BK7		0
I.1391.2541	meniscus	25,4		79,835/- 28,04		BK7	AR @ 1064 nm	3
xxxxxxxxx	meniscus	26			3,1			3
I.0111.3791	meniscus	37,9		109,5/-69,0	1,5	B270	AR @ 450-650 nm	2
I.0111.4002	meniscus	40	118,02	44,770/- 160,000	6,5	BK7		1
I.0111.4201	meniscus	42	56,03	131,520/- 36,310	3,5	BK7		0
I.0111.4321	meniscus	43,2	?	60,12/- 100,93	4,5	BK7	AR @ 1550 nm	1
I.0111.4322	meniscus	43,2	139,32	- 0,6934749 62	5	N-SF8	AR @ 1550 nm	3
xxxxxxxxx	meniscus	50		230,70/- 51,221	4	SF10		2
I.0110.6002	meniscus	60		- 0,8623188 41		BK7		4
xxxxxxxxx	meniscus	64,5	183	50,89/- 99,34	11			10
I.0111.8201	meniscus	82	196,29	63,310/- 158,000	11,5	BK7		0